

# Cardiac organoid cryopreservation kit

## Cardiac organoid cryopreservation kit

Cat. No. : RIPO-HWM005

### Description

The cardiac organoid cryopreservation kit is an optimized cryopreservation solution that allows whole cardiac organoid freezing at -80 °C and fast recovery. The kit supports high freeze/thaw viabilities with maximize function preserved. One kit allows the cryopreservation of 100 organoids.

### Product Specification

This cryopreservation kit is a serum-free, well-defined medium set with minimal batch variation. The cryopreservation medium contains DMSO.

### Product Information

Name	Size	Shipment	Storage
Cardiac organoid cryopreservation medium	10ml*3	Dry ice	The unopened products are stable for 6 months from the date of manufacture if stored under -80°C. The thawed products can be stored under 2-8 °C for two weeks.
Cardiac organoid recovery medium	10ml*3	Dry ice	

### Materials Required for Organoid Culture

- Ultra-Low Adherent 6 Well plate
- Programmed cell freezing container
- Cell freezing tube
- Human iPSC-Derived Cardiac Organoid Maintenance Kit (Cat. RIPO-HWM004)

### Equipment Required

- Incubator (37°C, 5% CO<sub>2</sub>)
- Orbital shaker (2 mm shaking diameter)
- Biosafety cabinet

- Refrigerator -80 °C

## Instruction of Use

### Cryopreservation

**Note: the cryopreservation can only be realized on cardiac organoids that have already cultured to the Maintenance stage (11 days+ of Acrobiosystems culture system).**

- Transfer each cardiac organoid into a cell freezing tube. Aspirate out the culture medium transferred with the organoid as possible and avoid causing damage to the organoid.
- Add 200ul of cryopreservation medium into each tube.
- Put the tubes into a programmed cell freezing container.
- Put the programmed cell freezing container into the refrigerator at -80 °C, keep overnight.
- Transfer the tube to other containers if necessary.

Note: The cryopreserved organoids can be stored under -80 °C for one year. Liquid nitrogen storage is not recommended.

### Recovery

- Take the cryopreserved organoids out from the refrigerator. Put the tubes into a water bath of 37°C for 1 min. (Shake the tubes gently with hands to accelerate thawing the organoids.)
- Verify the organoids are fully thawed.
- Transfer the organoids into the Ultra-Low Adherent 6 Well plate with in maximum 24 organoids per well.
- Add 5 ml of recovery medium in each well.
- Put the plate on an orbital shaker (as shown figures) with the speed of 100 rpm. Incubate at 37°C, 5% CO<sub>2</sub> for 48h.



### Culture

- After 48 h of recovery, change the recovery medium in each well to 5ml Cardiac organoid maintenance medium (Cat. RIPO-HWM004) per well

- b. Keep the plate on an orbital shaker (as shown in figures) with the speed of 100 rpm. Incubate at 37°C, 5% CO<sub>2</sub>.
- c. Full medium change every other days.

**Note: Organoids cannot be passaged.**

### Related Products

Product	Cat. No.
Human iPSC-Derived Cardiac Organoid Maintenance Kit	RIPO-HWM004