Product Data Sheet (DS)





NeuroFluidics NeoBento Dualink SHIFT LIGHT (Acro Certified)

Catalog No.: NFDLS-1

NeuroFluidics Line

High-throughput compartmentalized organs-on-chip devices for 2D cell culture

- Variety of architectures and related applications
- Variety of readouts compatibility
- Microchannels compartmentalization

Features

Specially designed to recreate synapses isolation and monitor axonal growth kinetics.

- LIGHT Version: 8 Chips & Data points per plate
- Discontinuous connectivity

Technical Specifications

• Asymmetrical shape with microchannels of different length (shorter length between the first and second channel than between the second and third channel).

reciment operations	
Surface Area:	• Channel 1: 18800 × 1000 × 200 μm (L × W × H), 18.8 mm² (32.9 mm² with reservoirs)
	• Channel 2: $6000 \times 200 \times 200 \mu m$ (L × W × H), 1.2 mm ² (15.3 mm ² with reservoirs)
	• Channel 3: 18800 × 1000 × 200 μm (L × W × H), 18.8 mm² (32.9 mm² with reservoirs)
	• Microchannels Tunnels: $100 \times 6 \ (\pm 1) \times 3,2 \ \mu m \ (L \times W \times H)$ for channel 1 to 2; $500 \times 6 \ (\pm 1) \times 100 \times$
	$3.2 \mu m (L \times W \times H)$ for channel 2 to 3; n=200; spaced by 20 μm
Volumes:	• Channel 1: 3.8 μL (117.7 μL with reservoirs)
	• Channel 2: 0.24 μL (114.1 μL with reservoirs)
	• Channel 3: 3.8 μL (117.7 μL with reservoirs)
Materials:	Microfluidic chip: PolyDiMethylSiloxane biocompatible and low compound absorbing (layer)
	170 μm thick + refractive index: 1.4)
	• NeoBento: Polystyrene (1.4 mm thick + refractive index: 1.59)
Formats:	• Microfluidic chip: 3 × 2 wells
	• QuarterBentos: 4 chips $(52,6 \times 34,6 \times 6,2)$
	• NeoBento: SLAS standard 96-well plate (127,8 \times 85,5 \times 17,1 mm)
Functions and R	deadouts deadouts
Capabilities:	Co-culture & compartmentalization
	hiPSC derived cell
	Synaptic isolation
	Functional analysis
Applications:	Study of synapses (pre-, post- and synaptic compartment)
	Synaptic transmission and localization
	Mitochondrial transport
	Microglial cells migration
	Neuroinflammation
	• Immunofluorescence
Doodouts	Immunofluorescence
Readouts:	Immunofluorescence Live Dead Assays
Readouts:	Immunofluorescence Live Dead Assays Live Staining



Product Data Sheet (DS)





• ELISA
• Calcium Imaging
• Electrophysiology

Acro Certify Disclaimer

This product is one of ACROBiosystems' Acro Certify products. ACROBiosystems and our Acro Certify partners have established a close partnership that includes an in-depth review of quality management and quality audits this product. Products from our Acro Certify partners have been qualified by ACROBiosystems to be included under Acro Certify. ACROBiosystems may provide Product information, including technical information, specifications, recommendations, literature, and other material (collectively, "Product Information") for customer's convenience. The accuracy and completeness of Product Information is not guaranteed and is subject to change without notice. ACROBiosystems is not responsible for the intellectual property or impact to intellectual property for products sold under Acro Certify.