

# CD Feed X Supplement Plus, Liquid

**Cat. # CF-1116-12**

**For Research Use Only**

**CD Feed X Supplement** is chemically defined, dramatically boost cell density and prolong cell viability thus increase expression of target protein, 5-10%(V/V)dosage in use. Used with CD 293 TGE Medium, BPM SBPfectin and 293 Expression MAX™-1 together is recommended for higher yield of protein expression.

## Product Features:

Significantly increase the productivity  
 Animal Origin-Free, Protein-Free, Serum-Free and Chemically Defined

## Storage:

Store **CD Feed X Supplement** at 2°C - 8°C in the dark.

## General Specification:

Category	Parameters
Form	Liquid
Volumn	500mL
pH	7.2 to 7.4
Endotoxin Testing	<10 EU/mL
Osmolality	950 to 1100 mOsm/kg
Mycoplasma Testing	Negative
Fungus Testing	Negative
Bacterial Testing	Negative

For detailed information or for further support, please contact our technical support team at:  
**[TechSupport@acrobiosystems.com](mailto:TechSupport@acrobiosystems.com)**