



Synonym

5-HT2A, HTR2A, HTR2, 5-hydroxytryptamine receptor 2A, 5-HT-2, 5-HT-2A, Serotonin receptor 2A

Source

Human 5-HT2AR Protein (S162K, M164W, S372N), Flag,His Tag(5HR-H55D3) is expressed from Baculovirus-Insect cells. It contains AA Asp 2- Ala 265 & Gln 313 - Lys 403 (Accession # [P28223-1](#) (S162K, M164W, S372N)). Predicted N-terminus: Asp

Molecular Characterization

This protein carries flag tag at the N-terminus and polyhistidine tag at the C-terminus.

The protein has a calculated MW of 55.3 kDa.

Endotoxin

Less than 1.0 EU per μg by the LAL method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Supplied as 0.2 μm filtered solution in 50 mM HEPES, 150 mM NaCl, pH7.5, CHS with glycerol as protectant.

Contact us for customized product form or formulation.

Shipping

This product is supplied and shipped with dry ice, please inquire the shipping cost.

Storage


Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- The product MUST be stored at -70°C or lower upon receipt;
- -70°C for 3 months under sterile conditions.

Bioactivity-ELISA

Human 5-HT2AR Protein (S162K, M164W, S372N), Flag,His Tag ELISA
0.5 μg of Human 5-HT2AR Protein (S162K, M164W, S372N), Flag,His Tag per well



Immobilized Human 5-HT2AR Protein (S162K, M164W, S372N), Flag,His Tag (Cat. No. 5HR-H55D3) at 5 $\mu\text{g/mL}$ (100 μL /well) can bind Anti-5-HT2AR Antibody (Cat. No. 5HR-H55D3) with a linear range of 0.02-0.625 $\mu\text{g/mL}$ (QC tested).

Background

This gene encodes one of the receptors for serotonin, a neurotransmitter with many roles. Mutations in this gene are associated with susceptibility to schizophrenia and obsessive-compulsive disorder, and are also associated with response to the antidepressants.

Clinical and Translational Updates

Discounts, Gifts,
and more!



Human 5-HT2AR Protein (S162K, M164W, S372N), Flag,His Tag (Detergent)

Catalog # 5HR-H55D3



Discounts, Gifts,
and more!



www.acrobiosystems.com

6/27/2024