

Biotinylated Monoclonal Anti-Mouse IL-4 Antibody, Rat IgG1 (AS542)

Catalog # FABr015-B



Application

Flow Cytometry.

Specificity

This product is a specific antibody specifically reacts with mouse IL-4 protein.

Immunogen

Purified mouse IL-4 Protein.

Conjugate

Biotin

Formulation

Supplied as 0.2 μ m filtered solution in PBS, pH7.4 with BSA as protectant.

Contact us for customized product form or formulation.

Storage

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- Store at 2° to 8°C and protect from direct light.
- Do not use after the expiration date shown on the label.

Background

Detection antibody for Multiplex Bead Assay.

Clinical and Translational Updates

Discounts, Gifts,
and more!

